

## CE-TEK 4000 All Weather Outdoor AED Enclosure

- IP66 standard, sealed against high pressure water hoses and dust.
- Fully corrosion resistant in UV stabilised high impact polycarbonate - the same material used in Police riot shields.
- Marine grade stainless steel keypad lock, salt spray tested.
- Suitable for coastal locations.
- Internal LED motion sensing light which can be removed and used as a torch during deployment.
- Thermostatically controlled heater ensures the defibrillator is kept in optimum conditions.
- Internal hook accommodates any FDA approved defibrillator.
- Large viewing window allows easy checking on the status of the AED defibrillator.
- Internal hinges that do not allow for external tampering.
- Double skin clamshell design provides additional insulation in cold weathers.
- Comes complete with comprehensive installation instructions and a fixing kit.
- Cabinet dimensions 18.3"x17.1"x9.4".
- **Manufactured in the UK**



**10 year warranty (5 year warranty on heater and keypad lock)**